

# **SOP on Histology Slide Preparation**

**Department of Anatomy**

**Faculty of Medicine**

**UWUSL**

# **Standard Operating Procedure (SOP) on Histology Slide Preparation**

**Issued By,**

Department of Anatomy, Faculty of Medicine, Uva Wellassa University of Sri Lanka.

## **Purpose**

To ensure uniform preparation of high-quality histology slides for teaching and laboratory purposes.

### **(1) Scope**

Applicable to histology laboratories and pathology laboratories.

### **(2) Materials Required**

#### **(A) Chemicals**

- 10% formalin (neutral buffered)



- Ethanol solutions



- Xylene



- DPX mounting medium



- Paraffin wax




- Hematoxylin and eosin stains (H & E)



## **(B) Equipment**

- Tissue processor



- Rotary microtome



- Tissue embedding machine



- Water bath



- Glass slides



- Cover slips



- Plastic cassettes



- Scalpels



- Plastic pipettes



- Labels and markers



### (3) Procedure


#### Step 1: Fixation – Preserve tissue in 10% formalin.

- The tissue specimen is immediately placed in 10% formalin to preserve its structure and prevent decomposition.
- Fixation stabilizes cellular proteins and prevents autolysis and bacterial growth.
- Proper fixation is essential for maintaining tissue morphology during further processing.




#### Step 2: Dehydration & Clearing

- The fixed tissue is passed through ascending concentrations of alcohol using a tissue processing machine, usually from 70% to 100% ethanol.
- This process removes water from the tissue.
- Gradual dehydration minimizes tissue shrinkage and distortion.



- After dehydration, alcohol within the tissue is replaced with xylene or another clearing agent using the same tissue processing machine.
- Xylene makes the tissue transparent and prepares it for wax infiltration.
- This step is necessary because paraffin wax is soluble in xylene but not in alcohol.



**Step 4: Embedding – Embed tissue in paraffin wax using the tissue embedding machine.**




- The cleared tissue is infiltrated with molten paraffin wax and then embedded in a wax block.
- Embedding provides support and firmness to the tissue for thin-section cutting.
- Proper orientation of the tissue during embedding is important for accurate microscopic examination.



**Step 5: Sectioning – Cut thin sections using a microtome.**

- The paraffin block is trimmed, and thin tissue sections are cut using a microtome, usually at 4–5  $\mu\text{m}$  thickness.
- Thin, even sections are essential for clear microscopic visualization.



### Step 6: Immerse in the water bath

- After sectioning with a microtome, the paraffin ribbon is gently placed on a warm water bath maintained at about 40–45°C.
- The warmth helps the tissue sections flatten and removes wrinkles or folds.
- Once the section spreads properly, it is carefully picked up onto a clean glass slide.
- The slide is then dried on a hot plate or in an incubator at 37–60°C for about 15–30 minutes to ensure proper adhesion and to remove excess water.



### Step 7: Staining – Perform hematoxylin and eosin staining.

#### (A) Deparaffinization

Sections mounted on glass slides are placed in xylene to completely remove paraffin wax. This step makes the tissue accessible to aqueous stains.

Time:

- Xylene I – 5 minutes
- Xylene II – 5 minutes



#### (B) Hydration

The slides are passed through descending grades of ethanol to rehydrate the tissue gradually. Finally, the slides are rinsed in water.

Time:

- 100% ethanol – 2 minutes
- 95% ethanol – 2 minutes
- 70% ethanol – 2 minutes
- Running tap water – 2 minutes



### (C) Hematoxylin Staining

Slides are immersed in hematoxylin.

Hematoxylin is a basic stain that binds to acidic structures such as the nucleus and nucleic acids, staining them blue-purple. This helps identify nuclear details clearly.

Time:

- Hematoxylin stain – 5–10 minutes



### (D) Washing and Bluing

Excess hematoxylin is washed off with water. Bluing agents or alkaline tap water convert the reddish-purple hematoxylin to a crisp blue color.

Time:

- Tap water wash – 2 minutes
- Bluing step – 1 minute



### (E) Eosin Staining

Slides are stained with eosin, which colors the cytoplasm, connective tissue, and extracellular components pink.

Time:

- Eosin stain – 1–2 minutes



### **(F) Dehydration**

The stained slides are dehydrated again using ascending alcohol concentrations to remove water from the tissue.

Time:

- 70% ethanol – 1 minute
- 95% ethanol – 1 minute
- 100% ethanol – 2 changes, 1 minute each



### **(G) Clearing**

Alcohol is replaced with xylene, making the tissue transparent and ready for mounting.

Time:

- Xylene I – 2 minutes
- Xylene II – 2 minutes




### **(H) Mounting**

A drop of DPX mounting medium is placed on the section and covered with a cover slip. The slide is then allowed to dry.

Time:

- Mounting procedure – 2–3 minutes
- Drying time – several hours to overnight



### Approximate Total Time


The complete H&E staining procedure usually takes about 30–45 minutes, excluding drying time.

- Staining enhances contrast and allows different tissue structures to be identified under the microscope.



### Step 8: Mounting – Apply coverslip with mounting medium.

- After staining, a coverslip is placed over the tissue section using DPX mounting medium.
- Mounting protects the stained section from damage and preserves it for long-term storage.
- Proper mounting also improves optical clarity during microscopic examination.



### (4) Quality Control

- Ensure tissue sections are wrinkle-free.
- Confirm adequate staining quality.
- Label slides correctly.

## **(5) Safety Precautions**

- Wear gloves and a lab coat.
- Handle xylene and formalin in a ventilated area.
- Dispose of sharps safely.

## **References**

1. Mescher AL. *Junqueira's Basic Histology: Text and Atlas*. 17th ed. New York: McGraw-Hill Education; 2024.
2. Young B, O'Dowd G, Woodford P. *Wheater's Functional Histology: A Text and Colour Atlas*. 7th ed. Philadelphia: Elsevier; 2023.
3. Kiernan JA. *Histological and Histochemical Methods: Theory and Practice*. 5th ed. Banbury: Scion Publishing; 2021.

### **Edited by,**

#### **Dr. Lakmal Hewage**

MBBS (COL), MD SURGERY (COL), MRCS (ENG)  
Consultant Neurosurgeon, Senior Lecturer (Grade II)  
Head of the Department of Anatomy

#### **Snr. Prof. Muditha Vidanapathirana**

MBBS (Col), DLM, MD, MA (SJP), FFFLM (UK)  
Dean, Faculty of Medicine, Uva Wellassa University of Sri Lanka

#### **Mrs. Nishadi Sachinthani Rodrigo**

BSc. Anatomy (Special) (USJ), MSc (UOC), PG. Cert. (Andrology) (UOC)  
Lecturer (Probationary)

#### **Dr. H.M. Piyumika Prabodhani**

MBBS (Ruhuna)  
Temporary Demonstrator, Department of Anatomy

#### **Dr. M.M. Chathurika Piyumali**

MBBS (Peradeniya)  
Temporary Demonstrator, Department of Anatomy

#### **Mr R. Kulendran**

Chief Technical Officer  
Department of Anatomy, Uva Wellassa University of Sri Lanka.