

SOP

**Acid Fast Staining Using Ziehl–Neelsen
Technique**

Department of Microbiology

Faculty of Medicine

UWUSL

Standard Operating Procedure (SOP)

Acid Fast Staining Using Ziehl–Neelsen Technique

Title:

Identification of Acid-Fast Bacilli Using the Ziehl–Neelsen Staining Technique

Issued By:

Department of Microbiology
Faculty of Medicine, Uva Wellassa University

(1) Purpose

To establish a standardized procedure for performing Ziehl–Neelsen staining in microbiology practical sessions to identify acid-fast bacilli (AFB), particularly *Mycobacterium* species, based on their ability to retain carbol fuchsin after acid-alcohol decolorization.

(2) Scope

This SOP applies to all medical students, academic staff, and laboratory personnel involved in microbiology practical sessions, including:

- Preparation of bacterial smears
- Performing Ziehl–Neelsen staining
- Microscopic identification of acid-fast bacilli

(3) Responsibilities

(3.1) Students

- Follow correct staining procedures
- Handle reagents and equipment safely
- Accurately observe and record findings

(3.2) Demonstrators / Lecturers

- Supervise staining procedure
- Ensure correct technique is followed
- Assist in interpretation of results

(3.3) Laboratory Staff

- Provide reagents and materials
- Ensure proper functioning of equipment
- Maintain laboratory safety standards

(4) Principle

- Acid-fast bacteria contain a high concentration of mycolic acids in their cell walls, making them impenetrable by ordinary aqueous-based staining solutions and difficult to stain by Gram's method, although they are structurally Gram-positive.
- The lipid-rich cell wall of acid-fast organisms stains with carbol fuchsin and resists decolorization with acid-alcohol. Therefore, these organisms retain the red/pink dye even after treatment with acids and alcohol.
- Only Mycobacteria with their thick waxy coats resist decolorization and remain stained with carbol fuchsin. These organisms are referred to as "acid-fast bacilli (AFB)" or "acid-alcohol-fast bacilli."
- The smear is counterstained with methylene blue, which stains the background and non-acid-fast organisms blue.

(5) Equipment and Materials

- Sterile inoculating loop
- Strong carbol fuchsin
- Acid-alcohol decolorizer
- Methylene blue
- Distilled water
- Bunsen burner
- Glass slides
- Staining rack
- Compound microscope
- Immersion oil

(6) Procedure

Step 1: Smear Preparation & Heat Fixing

Label one end of a clean, dry glass slide using a grease pencil to ensure proper identification. Prepare a thin smear of the specimen on the slide and allow it to air dry completely.




Pass the slide briefly through a flame 2–3 times to heat-fix.

Explanation: Heat fixing kills microorganisms, helps them adhere to the slide, and preserves cellular morphology during staining.

Step 2: Application of Strong Carbol Fuchsin

- Flood the smear completely with strong carbol fuchsin




- Heat intermittently to keep the slide steaming for approximately 5 minutes. Do not boil. Add more stain if necessary to prevent the slide from drying.
- Wash the slide gently with water.




- Explanation: Heating facilitates penetration of carbol fuchsin through the thick waxy lipid-rich cell wall containing mycolic acids.

Step 3: Decolorization with Acid-Alcohol

- Wash the slide with acid-alcohol until the color disappears from the smear.




- Immediately rinse with water.



- Explanation: Acid-alcohol removes the stain from non-acid-fast organisms, while acid-fast organisms resist decolorization and retain the red/pink stain.

Step 4: Counterstaining with Methylene Blue

- Flood the smear with methylene blue for 30 seconds.




- Wash with water and allow the slide to air dry.



- Explanation: Methylene blue stains non–acid-fast organisms and the background blue, providing contrast for easier visualization of acid-fast bacilli.

Step 5: Microscopic Observation


- Place the stained slide on the microscope stage.
- Focus under low power (10x) to locate the smear.
- Switch to high power (40x).
- Add a drop of immersion oil and examine under oil immersion (100x).



Observe and record:


- Presence or absence of AFB
- Morphology and staining characteristics

(7) Results Interpretation



Morphology of AFB

- Pink bright red bacilli
- Slightly curved appearance
- Beaded appearance
- Blue background



(8) Precautions

- Use a thin smear for proper staining
- Do not boil during heating
- Prevent stain from drying during steaming
- Do not over-decolorize
- Use clean slides and fresh reagents
- Follow aseptic technique
- Handle flame carefully

(9) Post Practical Procedure

- Dispose of slides safely
- Clean microscope stage
- Turn off light source
- Clean work area
- Dispose of contaminated materials according to laboratory safety guidelines

(10) Records and Documentation

- Practical number and title
- Date of experiment
- Observations and interpretation
- Student name and signature
- Demonstrator verification

(11) Reference

Source: National Institutes of Health (.gov) <https://share.google/2jy5Fi3oxBEIUk64q>

Prepared by:

Dr. Subodha Wickramasinghe
Head of the Department of Microbiology
Faculty of Medicine
University of Ruhuna

Dr. Lakmal Hewage
Acting Head of the Department of Microbiology
Uva Wellassa University of Sri Lanka

Snr. Prof. Muditha Vidanapathirana
Dean
Faculty of Medicine
Uva Wellassa University of Sri Lanka

Dr. Nadeesha Senevirathne
Temporary Demonstrator
Department of Microbiology

Submitted by:

Dr. Nadeesha Senevirathne
Temporary Demonstrator
Department of Microbiology

Date of Implementation: 06/05/2026